



“We were able to obtain and amplify DNA right there in the field. [...] miniPCR is ideal for the field.”

Elaine Guevara, Ph.D. Candidate at Yale University

miniPCR™: Studying genomes in the Madagascar forest

Whether due to bureaucratic red tape or the remoteness of the site, field researchers have long been plagued by the problem of how to bring the samples they collect back to a place where they can study them. Elaine Guevara was no exception to this rule. As a PhD candidate in biological anthropology at Yale, she made regular trips to the dry, spiny forests of Madagascar to visit a population of wild lemurs that has been the subject of long-term research by multiple groups, including her own. In particular, Guevara focused on the population’s genetic history and diversity, which she studied by capturing the animals and taking cheek swabs once in their lifetimes—a frequency (or, rather, infrequency) that made each sample precious.

But the colony of sifaka lemurs Guevara studied at Beza Mahafaly Reserve in Southwest Madagascar is endangered, and the logistical gymnastics required to transport biological matter belonging to such a species were complex. By the time researchers were able to remove the samples from the distant island nation and send them to a lab where Guevara would be able to extract and amplify the DNA, the swabs had already spent long periods stored under less-than-optimal conditions, jeopardizing both their own quality and the analysis that could be performed on them.



Faced with such an obstacle, Guevara considered an unusual alternative. If the difficulty lay in bringing the samples to the lab, she asked, why not bring the lab to the samples? She had heard from a colleague about a portable thermal cycler produced by the biotech startup Amplyus, and it seemed promising as a solution to her problem—a machine that she could carry into the field to amplify DNA on the spot. Powering the miniPCR through a lithium ion battery that she charged from a small solar panel, Guevara was able to begin performing whole genome amplification (WGA) assays straight from the remote Madagascar field site.

Only one question now remained: would the WGAs work? That is, would the miniPCR be able to stand up to the task at hand as well as the bench-top machines upon which the team had previously relied? The answer, Guevara asserts, is yes. “The miniPCR seems just as effective as a traditional thermocycler, but of course far more transportable, durable, and energy-efficient,” she says. It would seem that with the miniPCR, Guevara had found a way to overcome the challenge that had for so long stymied her peers—and the resulting enhancements in the condition of the DNA and the extent of analysis that could be conducted on it, of course, speak for themselves.

Today, having returned from Madagascar after a first round of gene analysis, Guevara is excited for the horizons that have opened to her research, and she hopes to continue to use the miniPCR to facilitate her

To learn more please visit our website www.minipcr.com or email us: team@minipcr.com

© 2016 by Amplyus, Cambridge, MA, USA. All rights reserved.



findings. “The use of the miniPCR represented a major improvement for us because we were able to obtain and amplify DNA right there in the field,” she says. “Going forward, I am hoping the miniPCR will also allow us to do even more of the genetic analyses in Madagascar.... I think miniPCR is ideal for the field.”

Elaine’s microsatellite genotyping of these lemurs will help determine relatedness for the population, and can answer questions about mating strategies, fitness, and dispersal, as well as help reconstruct the demographic history of the population and assess their genetic diversity.

Victoria Lin

About miniPCR™

miniPCR™ is a groundbreaking miniaturized thermal cycler offering fully-featured PCR performance for a fraction of the price.

- **Powerful:** miniPCR thermal cyclers deliver the same high-end performance as benchtop PCR platforms.
- **Portable:** miniPCR can fit in the palm of a hand and communicates seamlessly with your laptop, smartphone or tablet. Portable battery packs allow for on-site application in remote locations.
- **Simple:** Our unique, easy to use, software interface allows users to control and visualize their PCR reactions from their mobile device.
- **Open:** miniPCR is fully compatible with standard PCR tubes and reagents, making it ready to run existing assays.

miniPCR is based at the Harvard Launch Lab, an incubator for innovative companies run by Harvard University alumni based in Boston, Massachusetts, USA.

Feature	Specification
Sample format	8 x 0.2ml PCR tubes (strip compatible)
Max heating ramp rate	3.2°C / sec
Heated lid	Independent lid heater up to 120°C with PID control
Max cooling ramp rate	2.2°C / sec
Temperature control	Resistive heating; forced air cooling
Control system	Embedded thermistors and PID algorithm
Dimensions	2" x 5" x 4"
Weight	1 lb. (0.450 gr)
Power supply	AC 100-240V, 50-60hz, 70W
Battery operation	4h – 6h uninterrupted run time on Li-Ion battery

The miniPCR team of molecular biologists, engineers, and designers is dedicated to miniPCR’s mission of making science accessible to everyone, everywhere.

To learn more please visit our website www.minipcr.com or email us: team@minipcr.com

© 2016 by Amplyus, Cambridge, MA, USA. All rights reserved.