



P51™ Chlorophyll Lab: Green Glows Red!



For use with P51™ Molecular Viewer



(or other blue light illuminator¹)

Instructor's Guide Contents

- | | |
|--|-------|
| 1. Background and significance | p. 2 |
| 2. Laboratory guide | p. 4 |
| 3. Study questions | p. 10 |

¹ Compatible with blue light transilluminators such as blueGel™, blueBox™ and other 460-480 nm illuminators.

1. Background and significance

When we start a fire for warmth and burn wood, we are releasing energy stored in the bonds of plant molecules. When we look for nourishment by eating sugar or other carbohydrates, we are doing the same thing, releasing stored chemical energy from a plant for our own use. But energy can never be created nor destroyed, meaning the plant that we are using for energy had to obtain that same energy from somewhere else. So where did the plant matter we use for heat and nourishment get its energy in the first place? Plants obtain energy from the sun, and store that energy in the bonds of carbohydrates. In fact, virtually all energy used by any living organism on Earth started out as solar energy captured by plants or other organisms that perform photosynthesis. This ability to convert light energy into usable chemical energy provides nearly all of the energy used by life on this planet.


Photosynthesis is the process whereby light energy is absorbed and converted into chemical energy that can be used by organisms. The key first step of photosynthesis, therefore, is to efficiently absorb light energy. In plants and algae, this job is done primarily by the molecule *chlorophyll*. Chlorophyll is a pigment, a molecule that absorbs light. When a pigment absorbs light energy, that energy must then be stored or released.

When light hits a pigment in photosynthesis, the energy from the light is absorbed, and it causes an electron of the pigment molecule to raise in energy level. But that electron can't stay in that excited state forever; it will quickly pass that energy off, returning to its original state. If you know what it feels like to wear a black shirt on a bright sunny day, you know one place that energy can go; the dark pigments of the black shirt absorb the light energy but then release it as heat. The plant pigment chlorophyll will absorb light energy like any pigment, but instead of releasing that energy as heat, it will pass that energy from light off to other components of the chloroplast, where it will eventually be stored as chemical energy.

Chlorophyll located in the chloroplast thylakoid membrane is arranged in such a way that, when light is absorbed, the energy in the electron will be passed down a chain from one chlorophyll molecule to the next until it reaches the photosystem reaction center. In the reaction center, the energy is transferred down the electron transport chain, eventually contributing to the production of sugar molecules in the cell. In this way, energy enters the plant leaf in the form of light, but is transformed into stored chemical energy.

There are two main forms of chlorophyll, *chlorophyll a* and *chlorophyll b*. Both chlorophylls best absorb wavelengths in the blue to violet ranges of the spectrum and also in the red ranges of the spectrum, but each is best at absorbing slightly different wavelengths of light. Chlorophyll's inability to absorb green wavelengths is the reason plants are green—when white light from the sun hits a plant, the short violet and blue wavelengths are absorbed along with the long red wavelengths, but the intermediate green wavelengths are reflected. Plants also have other pigments that help absorb some of this light that chlorophyll misses. Carotenoids are yellow to red pigments present in almost every leaf that absorb some of the green that chlorophyll can't. It is these carotenoids that give fall foliage much of their yellow and orange colors.

Light energy that is absorbed by a pigment can be released as heat energy, as in a black shirt on a sunny day, or it can be transferred into chemical energy as in chlorophyll in a leaf, but there is another way that this energy can also be released. In some cases, absorbed light energy can be released back again as light. When this happens, because some energy



has been lost to heat and other processes, the light that is released is of a longer, lower-energy wavelength than when it was first absorbed. This is known as the Stokes shift and it is the basis of the phenomenon known as fluorescence. Chlorophyll absorbs light best in the violet to blue areas of the visible spectrum. This blue and violet light will be absorbed by chlorophyll, but if the chlorophyll cannot transfer the absorbed energy to be stored as chemical energy, that energy must go somewhere: it is released back as light of a longer wavelength—red light.

So why aren't most plants emitting red light? In a living plant, fluorescing chlorophyll would be a problem, as energy that is released as light cannot be used to make sugar. A fluorescent red plant would be very bad at storing energy in the form of sugar and would essentially starve to death. In fact, measuring the amount of red fluorescence produced by chlorophyll in a plant is actually a tool that scientists use to measure the efficiency of the chloroplast electron transport chain. When the electron transport chain is operating at high efficiency, very little fluorescence is produced because a high proportion of the energy absorbed by chlorophyll is being efficiently transmitted into chemical energy. If more fluorescence is produced, less of the absorbed sunlight is being transmitted into chemical energy.

In this lab, you will be able to observe and separate the pigments found in a leaf. You will then observe the appearance of those pigments—both in their normal physiological surrounding within a leaf, as well as when they are separated from the leaf and dissolved in solution.

2. Laboratory guide



Acetone and isopropyl alcohol are highly flammable. Keep away from heat sources or open flame. Protect eyes and skin. Gloves and protective eyewear should be worn during the entirety of this lab. Plant extract will stain clothing and other materials. Be careful while handling.

1. Extracting plant pigments

For extraction we recommend 100% acetone or at least 90% isopropyl alcohol

1. **Place a few leaves in a mortar.**
 - Any green leaf, from baby spinach to pine needles will work.
 - If using baby spinach 1-3 leaves should be plenty, depending on the size of your mortar.
2. **Add isopropyl alcohol or acetone as a solvent.**
 - Add just enough that the leaves are wet.
 - Adding too much liquid can make crushing with the pestle more difficult.
3. **Use the pestle to grind and crush the leaves.**
 - Grind the leaves using the pestle until they are thoroughly crushed and very few intact pieces of leaf are left.
4. **Add some more solvent if needed.**
 - Add more acetone or alcohol until there is just enough liquid to be able to pour.
 - Adding too much liquid will cause the extract to become too dilute.



If using a centrifuge:

5. **Add at least 200 μ l of the ground leaf to a 1.7 ml microcentrifuge tube.**
 - Use a plastic transfer pipette or carefully pour some of the liquid into the tube.
6. **Using a microcentrifuge, spin for at least 20 seconds at 10,000 rpm.**
 - Make sure that the centrifuge is balanced with equal volumes in all tubes.
7. **Carefully remove tube from centrifuge. The top liquid phase is your leaf extract.**
 - There should be a transparent green liquid at the top of the tube and an opaque pellet of plant matter at the bottom of the tube.
 - The transparent green liquid is your plant extract.
 - Handle the tube carefully so you don't disturb the pellet.



If using filter paper:


5. Fold a piece of filter paper and place it inside a funnel.
6. Place the funnel in a beaker, test tube, or other collection container.
7. Pour the entirety of the liquid in the mortar into the filter paper.
 - Wait for the liquid to filter through the paper (a few minutes).
8. The liquid that collects in your beaker is your leaf extract.

2. Paper chromatography

Your teacher may elect to use slightly different paper chromatography protocol.

Obtain a clean dry piece of filter paper and cut it into a long strip.

- The strip should be between 1 and 3 centimeters wide and at least as tall as the beaker you intend to perform chromatography in.
 - Cut the tip of the paper into the shape of a point.
2. Use pencil to draw a line on the paper about 2-3 cm from the tip.
 - The pencil line should be perpendicular to the length of the paper.
 3. Use a toothpick to draw a line of leaf extract on top of your pencil line.
 - Dip toothpick in the leaf extract and drag it across the pencil line.
 - Repeat until the entire line is wet.
 - Do not get too wet, or the extract will bleed away from the line.
 4. Let the leaf extract dry and reapply until there is a sufficient amount of leaf extract dried on the paper.
 - The solvent will air-dry quickly – usually in ~15 seconds.
 - When dry, use a toothpick to reapply plant extract on the same line as in step 3.
 - When finished the line should be a solid green color.
 - 5-8 applications are usually enough, depending on the concentration of your leaf extract.
 5. Add chromatography solvent your beaker.
 - Use 70% isopropyl alcohol so that the total amount of liquid is no more than 1 cm deep.
 6. Rest a glass stirring rod, pencil or other similar object across the top of the beaker.
 7. Hang your paper off of the stirring rod so that the tip of the paper is in the solvent, but the line of dried plant extract is above the solvent.
 - Fold the paper so that it can hang from the rod. Alternatively attach it to the rod using tape.
 - The point of the paper should be in the alcohol, but the dried plant extract should be above the solvent.
 8. Wait until you see significant separation of pigments (15-20 minutes).



3. Observing pigment fluorescence

1. Obtain a 200 μ l tube and label the tube "PE" for leaf extract.
2. Add up to 50 μ l leaf extract to your tube.
 - If using a 200 μ l pipette, transfer 50 μ l to the tube.
 - If using a 20 μ l pipette, transfer 20 μ l several times.
 - If using a plastic transfer pipette, fill the tube $\frac{1}{4}$ to $\frac{1}{2}$ full.
3. Prepare controls.
 - Add approximately the same volume of your solvent without any leaf extract added. Label the tube "S" for solvent.
 - Tear or cut a small piece of leaf and place it in a tube alone.
 - Optional: Add green food coloring to water and dilute it until is a similar color to your plant extract.
4. Place the tubes in P51 or other blue light illuminator and record your results.



Tube	Color of solution in ambient light	Color of solution in blue light illuminator	Notes
1. Plant extract			
2. Solvent only			
3. Leaf only			
4. Food coloring			

4. Isolating pigments for fluorescence observation

For isolation we recommend 100% acetone or at least 90% isopropyl alcohol.

1. Remove the chromatography paper from the beaker.
2. Use scissors to cut distinct bands of pigment out of the chromatography paper.
 - As best as you can, cut strips of paper so that each strip contains only a single color of extracted pigment.
 - The strips should be narrow enough to fit in a 200 μ l tube.
3. Place each strip of paper in a different 200 μ l tube.
 - Label the tubes 1, 2, 3 and so on for as many strips you have cut out.
 - #1 should correspond to the band on your paper that ran the farthest on your chromatography paper. #2 should have run the second highest and so on.
 - You will likely have about 2-3 discernible pigments, but as many as 5 are possible depending on your chromatography protocol.
4. Add ~50 μ l of acetone (or 90%+ isopropyl alcohol) to each tube.
5. Cap the tube and flick the tube several times to mix well.
 - The pigment should be washed from the paper and into the solvent in the tube.
6. Remove the paper from each tube.
 - You may use forceps, a toothpick or a pipette tip to pull the strip of paper out of the tube.
7. Place the tubes in P51™ or other blue light illuminator and record your results.



Tube#	Distance traveled on paper	Color of band on paper	Color of solution under blue light	Notes
1				
2				
3				
4				
5				

CER Table

Fill in the table based on your results from the lab.

Question:

Why do you observe red fluorescence under certain circumstances?

<p>Claim</p> <p>Make a clear statement that answers the above question.</p>	
<p>Evidence</p> <p>Provide data from the lab that supports your claim.</p>	
<p>Reasoning</p> <p>Explain clearly why the data you presented supports your claim. Include the underlying scientific principles that link your evidence to your claim.</p>	

	4	3	2	1
CLAIM <i>A statement that answers the original question/ problem.</i>	Makes a clear, accurate, and complete claim.	Makes an accurate and complete claim.	Makes an accurate but incomplete or vague claim.	Makes a claim that is inaccurate.
EVIDENCE <i>Data from the experiment that supports the claim.</i> <i>Data must be <u>relevant</u> and <u>sufficient</u> to support the claim.</i>	All of the evidence presented is highly relevant and clearly sufficient to support the claim.	Provides evidence that is relevant and sufficient to support the claim.	Provides relevant but insufficient evidence to support the claim. May include some non-relevant evidence.	Only provides evidence that does not support claim.
REASONING <i>Explain why your evidence supports your claim. This must include scientific principles/knowledge that you have about the topic to show why the data counts as evidence.</i>	Provides reasoning that clearly links the evidence to the claim. Relevant scientific principles are well integrated in the reasoning.	Provides reasoning that links the evidence to the claim. Relevant scientific principles are discussed.	Provides reasoning that links the evidence to the claim, but does not include relevant scientific principles or uses them incorrectly.	Provides reasoning that does not link the evidence to the claim. Does not include relevant scientific principles or uses them incorrectly.

Rubric Score	3	4	5	6	7	8	9	10	11	12
Equivalent Grade	55	60	65	70	75	80	85	90	95	100


We recommend that teachers use the following scale when assessing this assignment using the rubric.

Teachers should feel free to adjust this scale to their expectations.

3. Study Questions

Pre-lab questions

- The graph to the right shows the approximate relative absorbance of chlorophyll a and b. For reference, approximate color ranges of different wavelengths are listed below. If you were designing an extremely efficient light to grow plants, what wavelengths or colors would it emit? Justify your answer.



Violet: 380-450 nm	Blue: 450-495 nm
Green: 495-570 nm	Yellow: 570-590 nm
Orange: 590-620 nm	Red: 620-750 nm

- What wavelengths or colors would it be inefficient for your grow light to produce? Justify your answer.
- Plants also contain pigments known as carotenoids. Carotenoids also absorb light and pass that energy to chlorophyll. Carotenoids absorb light best in the range of about 450 nm to about 530 nm. Why might it be helpful for a plant doing photosynthesis to use carotenoids along with chlorophyll?
- Chlorophyll is a unique pigment in that it has evolved not just to efficiently absorb light, but also to transfer that energy to other molecules. Why is it so important that chlorophyll be good at transferring energy?
- Imagine that you left a leaf that was performing photosynthesis in the sun, and next to it, you left a material that was the same size, shape, mass, and color, but did not perform photosynthesis and contained pigments other than chlorophyll. And after half an hour you measured the temperature of both materials. Which material would you expect to be warmer? Justify your answer.

Post-lab questions

1. What color was the plant extract in ambient light, and what color was the extract in blue light?
2. For each control tube that used, explain why it was important to include that tube as part of your experiment. If you had not included that control in your experiment, what alternate conclusion could you have made?

Control	Without this control, what alternate conclusion regarding the red fluorescence could you make?
Solvent only	
Unextracted leaf	
Green food coloring	

3. Based on the graph of chlorophyll absorbance, is the red color you observed more likely from chlorophyll a or chlorophyll b? For reference, P51™ emits light with a wavelength of 465 nm. Explain your answer.
4. If you wanted to more clearly test the *other* chlorophyll molecule for fluorescence, what wavelength/color light would you use?
5. How many pigments did you see on your chromatography paper?
6. Of those pigments, which fluoresced and which didn't?
7. Red light emitted by the chlorophyll is a form of energy. Where did that energy come from?
8. If a pigment does not fluoresce, there are two possible reasons why. (1) It did not absorb the light, or (2) after absorbing the light it released the light energy in a form other than fluorescence. Can you propose a simple experiment to see if either of those two explanations is accurate?