


# GELATO™ electrophoresis and visualization system

QP-1600-01

## TECHNICAL SPECIFICATIONS



|                                |  |
|--------------------------------|--|
| <b>RUNNING VOLTAGE</b>         | ADJUSTABLE. 50, 75, 100 AND 135 V  |
| <b>ILLUMINATION</b>            | INTEGRATED, 120 HIGH-INTENSITY BLUE LEDS   |
| <b>BUFFER COMPATIBILITY</b>    | TAE AND TBE  |
| <b>GEL SIZE</b>                | SMALL TRAY 6 X 6 CM, LARGE TRAY 12 X 6 CM  |
| <b>LOADING CAPACITY</b>        | UP TO 26 LANES ON SMALL GELS AND 50 ON LARGE GELS                                |
| <b>MULTICHANNEL COMPATIBLE</b> | YES  |
| <b>COMPATIBLE DNA DYES</b>     | SEEGREEN™, GELGREEN®, SYBR™, GREENVIEW™ AND OTHER BLUE-LIGHT COMPATIBLE DNA DYES |
| <b>DOCUMENTATION</b>           | FOLD-A-VIEW™ PHOTO HOOD  |
| <b>TIMER</b>                   | VARIABLE, UP TO 120 MIN  |
| <b>DIMENSIONS</b>              | 25 X 17 X 8 CM   |
| <b>WEIGHT</b>                  | 800 GRAMS  |
| <b>POWER SUPPLY</b>            | DC 9 V, 4 A  |

© Ampliys LLC, all rights reserved

Designed in Cambridge, MA | Assembled in the USA | 1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com