


# Gyro™ Plus

Variable speed microcentrifuge  
QP-1800-02



## TECHNICAL SPECIFICATIONS

### FEATURES

Variable speed microcentrifuge with multiple speed settings and a user-friendly LED display and interchangeable rotors.

### SPEED RANGE

1,000 to 12,000 RPM

### MAXIMUM RCF

6,900 g

### ROTOR CAPACITY

6 x 1.5 ml, 16 x 0.2 ml, 6 adapters for 0.2 or 0.5 ml

### NOISE LEVELS

< 60 dB

### DIMENSIONS

17 x 15 x 13 cm

### WEIGHT

1 kg

### VOLTAGE

AC 85-240 V, 50/60 Hz

### ADDITIONAL FEATURES

LED display, automatic safety stop

1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com