


P51TM Molecular Fluorescence Viewer

QP-1900-01



TECHNICAL SPECIFICATIONS

ILLUMINATION	HIGH INTENSITY BLUE LEDs, ~470 nm WAVELENGTH
FILTERS	AMBER AND YELLOW, BLUE LIGHT FILTERING
POWER SOURCE	USB-C or 9V battery (NOT INCLUDED)
RACK CAPACITY	8 POSITIONS
TUBE FORMAT	0.2 ML TUBES, DOMED OR FLAT CAPS, STRIP COMPATIBLE
CONSTRUCTION	INNER ILLUMINATION CHAMBER, OUTER HOLDING CHAMBER
ASSEMBLED DIMENSIONS	120 MM X 90 MM X 60 MM
PACKAGE DIMENSIONS	240 MM X 33 MM X 200 MM
WEIGHT	153 G (WITH BATTERY)
COMPOSITION	50% POST-CONSUMER RECYCLED MATERIALS

Designed in Cambridge, MA | Assembled in the USA | 1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com