


Viewit™

QP-1401-01



TECHNICAL SPECIFICATIONS

FEATURES

- Transilluminator accessory for Bandit™
- Orange filter plate for viewing fluorescence
- High intensity blue LED strips
- Low voltage USB power

ILLUMINATION

2 strips with a total of 12 high-intensity blue LEDs at ~470 nm wavelength

STAIN COMPATIBILITY

SeeGreen™ - GelGreen® - GreenView™ Plus - EvaGreen - SYBR Safe - GRSafe, Gel Star - SYBR Gold, SYBR Green I and II - SYPRO Ruby and others compatible with blue-light illumination

INPUT VOLTAGE

5 V DC, USB-C

OUTPUT VOLTAGE

40 V DC

WEIGHT

3.56 oz (101 grams)

DIMENSIONS (L x W x H)

5.35 x 3.46 x 0.91 in (136 x 88 x 23 mm)

1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com