

blueBox™ Pro Transilluminator with Imaging Hood

QP-1700-03

TECHNICAL SPECIFICATIONS



ILLUMINATION

139 HIGH-INTENSITY BLUE LEDS AT ~470 nm WAVELENGTH

STAIN COMPATIBILITY

GelGreen™ - GreenView™ Plus - EvaGreen - SYBR Safe -

GRSafe, Gel Star - SYBR Gold, SYBR Green I and II - SYPRO Ruby

and others compatible with blue-light illuminator

VIEWING AREA

140 mm X 114 mm

CONSTRUCTION

100% SCRATCH-RESISTANT ACRYLIC

IMAGE CAPTURE

FOLD-A-VIEW™ PHOTO HOOD. LOCK-AND-CUT HINGE

FOR EASY BAND CUTTING AND GEL EXCISION.

INPUT VOLTAGE

100-240VAC, 50-60Hz.

OUTPUT VOLTAGE

12V DC

Designed in Cambridge, MA | Assembled in the USA | 1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com