

blueBox™ Pro Transilluminator with Imaging Hood

QP-1700-03

TECHNICAL SPECIFICATIONS



ILLUMINATION	139 HIGH-INTENSITY BLUE LEDS AT ~470 nm WAVELENGTH
STAIN COMPATIBILITY	GelGreen™ - GreenView™ Plus - EvaGreen - SYBR Safe - GRSafe, Gel Star - SYBR Gold, SYBR Green I and II - SYPRO Ruby and others compatible with blue-light illuminator
VIEWING AREA	140 mm X 114 mm
CONSTRUCTION	100% SCRATCH-RESISTANT ACRYLIC
IMAGE CAPTURE	FOLD-A-VIEW™ PHOTO HOOD. LOCK-AND-CUT HINGE FOR EASY BAND CUTTING AND GEL EXCISION.
INPUT VOLTAGE	100-240VAC, 50-60Hz.
OUTPUT VOLTAGE	12V DC
DIMENSIONS	6.125" (l) x 5.125" (w) x 4" (h) (15.5 cm x 13 cm x 10 cm)
WEIGHT	440 g

©2020 - minipcr bio, all rights reserved.

Designed in Cambridge, MA | Assembled in the USA | 1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com

For blueBox orders and purchasing options visit minipcr.com