

blueGel™

Integrated electrophoresis and visualization system
QP-1500-01



TECHNICAL SPECIFICATIONS

FEATURES

- Gel electrophoresis unit with built-in transilluminator
- Seamless polycarbonate and ABS injection molded parts
- Two-button operation; one for Run, one for Light
- Casting system with integrated comb holder
- Platinum and stainless steel electrodes
- High intensity blue led transilluminator panel

GEL SIZE

60 x 60 mm

GEL VOLUME

15 to 20 ml

BUFFER VOLUME

30 ml, optimized for TBE

ONE OR TWO ROWS OF COMBS

9 or 13 wells each

POWER SUPPLY

AC 100 - 240 V, 50-60 Hz, 48 V

DIMENSIONS

9" x 4" x 3" height (23 x 10 x 7 cm.)

WEIGHT

0.75 lbs., 350 gr

Designed in Cambridge, MA | Assembled in the USA | 1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com